

November 29, 2018

Ref: Approval for Use Statement: Nanosonics Trophon2 - Lumify

Philips has approved for use Nanosonics Trophon2 for the following Lumify transducer models: C5-2, L12-4, and S4-1, and future devices based on the Quick Disconnect USB architecture when used per label instructions.

NOTE: These transducers types, *Quick Disconnect USB architecture*, are not approved for use with the original Trophon device due to cable diameter limits.

This approval is based on both component level (historical testing) and finished device testing conducted at Nanosonics. Since it is impossible to predict with absolute certainty that these laboratory conditions mimic all field use conditions, approval for use can be withdrawn should field reliability of these devices become compromised.

Reference Nanosonics Test Reports
PJTMAT-REP 338 C5-2 (A493465)

NOTE: Updates to product literature on the Philips Transducer Care web site will occur at the next revision of these documents.

Kevin G. Wickline
Chemist - R&D Engineering
Philips Healthcare

